camerajet

powered by Gap Optical Targeting™

Instructions for use

dyson

Dyson CameraJet™

powered by Gap Optical Targeting™

Instructions for use

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


Descriptive Information

When the device should not be used

Where the following statements apply, please consult a dentist or a healthcare professional before using the device.

- Children below 18 years or as per local regulations.

Individuals with:

- Oral or facial implants or pacemakers.
- Severe periodontal disease or untreated dental infections.
- Diagnosed bleeding disorders or who are on anticoagulant therapy
- Mouth ulcers, open wounds, or severe stomatitis.
- Prior oral surgery, dental implants, or teeth extractions within the last 2 months.
- Oral or dental conditions that may be exacerbated by brushing or interdental cleaning.
- Reduced physical, sensory, or mental abilities, or those lacking experience or knowledge, should not use the device unless under supervision or having received instructions from someone responsible for their safety.

General warnings and precautions

Read this information carefully before use and save it for future reference.

- Only use this device as indicated in this booklet.
- **WARNING:** No modification of this equipment is allowed.
- If the device is damaged, stop using it immediately.
- Replace the brush head every 3 months, or sooner if damaged or excessively worn.
- Use only Dyson-manufactured brush heads, rechargeable batteries and water tanks.
- **WARNING:** Use only an external supply with the following specifications: Listed/certified USB-C power adapter with safety extra-low voltage (5V DC, 500mA–3A). Non-certified sources can be hazardous.
- Keep the charging contacts clean to effectively protect from oxidation and other possible damage.
- Only perform maintenance as described in this document or advised by Dyson.
- **WARNING:** Do not attempt to clean or repair the product while the device is active.
- Supervise children to prevent them from playing with the device.
- **WARNING:** Keep the charging cable away from toddlers and children to avoid potential hazard.
- For hygiene, each device is recommended for single person use only.
- **CAUTION:** Follow cleaning instructions in this booklet to prevent buildup of residue or biofilm in the fluid system.
- **⚠ CAUTION:** Do not direct the jet towards sensitive areas like the eyes, nose, or ears.
- Oral jewellery should be removed before using the device.
- Only use cool or warm liquid, Dyson mouth rinse where available, or an alternative low-foaming mouth rinse.
- Use only at altitudes up to 3000m.
- Do not expose the device to oils and solvents, including alcohol-based wipes, perfume, skincare products and essential oils.
- If excessive gum bleeding occurs or persists, discontinue use and consult a dentist.
- Avoid brushing immediately after consuming acidic food or drinks.
- Travel case shall not be used during usage of the device.

Battery information

The battery must be charged at an environmental temperature between 5 °C and 35 °C. Do not expose the device or battery to fire or temperatures above 70 °C.

This device has a sealed lithium-ion battery, which is safe under normal use. In the unlikely event that liquid leaks from the device, do not attempt to operate the device, do not touch the liquid as it may cause irritation or burns. If contact occurs, rinse the affected area thoroughly with water and seek medical attention immediately.

Maintenance and repair 1

It is recommended to replace your brush head every 3 months.

When replacing battery packs, use only the Dyson battery Part Number 975294-01 for this device. Other types of batteries are incompatible and may lead to severe failure.

To replace the battery **2**

1. Before removing the battery, disconnect the charging cable, and ensure the device is in standby by holding the Power button for 2 seconds. You can check the device is in standby by moving it. If any indication lights illuminate, the device is not in standby.
2. Remove the adhesive label from the base of the device, unscrew the screw using a T6 screwdriver, and slide out the battery from the brush handle.
3. Insert a new battery by aligning it with the corresponding U-shaped mark in the battery compartment at the base of the handle, so the connector fits with the correct polarity, and pushing the battery in.
4. Insert and tighten the screw, being careful not to overtighten, and paste the adhesive label onto the base of the battery.
 - **⚠ CAUTION:** Do not use a screw other than the screw provided.
 - **⚠ CAUTION:** When the battery is removed, keep battery terminals insulated and away from metal objects to avoid short-circuiting.
 - Dispose of the used battery in accordance with local regulations or guidelines.
 - Do not use modified or damaged batteries.
 - If the device is to be stored unused for a long period, the batteries should be removed.

All repairs other than those above can only be carried out by Dyson service personnel or a third party appointed by Dyson. Servicing and repair by any other party is prohibited.

Disposal

The battery should be removed from the device before disposal. See "Maintenance and repair" section above for removal instructions. Dispose of or recycle the battery in accordance with local regulations or guidance.

This symbol indicates that this device should not be disposed with your general household waste. To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources. When this device has reached the end of its life, please take it to a recognised WEEE (Waste Electrical and Electronic Equipment) collection facility such as your local civic amenity site for recycling. Your local authority or retailer will be able to advise you of your nearest recycling facility.

Operating instructions

Setup **3**

1. Remove the Travel case by pressing the release button at the bottom of the cover.
2. Remove the device from the dock and attach the brush head to the toothbrush handle.
Note: Both the brush head and handle should be facing in the same direction. Keep the camera dry before attaching the brush head for optimal performance.
3. To ensure optimal camera vision and gap detection keep the lens cover clean. If you experience limited jetting, pause the clean by pressing the Power button, remove the brush head, wipe the lens clean, replace brush head and press the Power button to resume clean.
4. Slide the tank up, without twisting, to expose the inlet.
5. Fill with Dyson mouth rinse where available, or an alternative low-foaming mouth rinse.
⚠ CAUTION: Do not jet with hot liquid.
6. Assemble the tank by pushing it down on to the brush handle.
7. If powered off, press the Power button and wait for 6-8 seconds to activate the camera. Upon first use, the device is automatically in **Brush and Auto-jet mode**. After this, it will remember your settings.
8. If in standby, press any button or move the device and wait 6-8 seconds for the camera to activate.

Modes

The device has three modes. If you wish to change modes, do so before activating the device.

Recommended: **Brush and Auto-jet mode** simultaneously brushes and Auto-jets. Auto-jet uses the camera to automatically detect interdental gaps and precision floss.

- **Auto-jet mode** auto-jets without brushing.
- **Brush mode** brushes without Auto-jet.

Each brushing mode has three intensities to choose from. Gentle, variable and deep clean. You can also trigger the jet at any time by pressing the Jet button.

Technique

Brush and Auto-jet mode:

1. Split the lower and upper jaw into 3 sections – left, right and centre.
2. Pick a section to clean. Start with the outer surface. Hold the toothbrush horizontally and move the brush head across the teeth in small, oval, back-and-forth motions.
3. Move onto the inner surface, and then the chewing surface, using the same brushing motion.
4. When brushing the inside of the front teeth, tilt the brush slightly upright and brush up and down on each tooth.
5. The pacer will vibrate the toothbrush to indicate when it is time to move on to the next section.

Auto-jet mode:

- Floss the outside and inside of the teeth.

Brush mode:

- Follow the same brushing technique described for Brush and Auto-jet mode.
-

Use 4

1. Facing the brush head into a sink, press and hold the Jet button until all air is expelled and the nozzle begins to jet fluid.
 2. If you are using Brush and Auto-jet mode, or Brush mode, apply toothpaste. Ensure the camera lens is not covered. If you are using Brush and Auto-jet mode, or Brush mode, apply toothpaste. Use the Dyson toothpaste where available, or a non-foaming alternative. Ensure the camera lens is not covered. Place the toothbrush head flat against your teeth. Hold horizontally to follow the jawline.
 3. To begin operation, press the Power button.
 4. Use light pressure on the teeth and gums and let the toothbrush do the work.
 5. Guide the toothbrush slowly across the teeth, moving from tooth to tooth with a gentle oval motion. Do not scrub.
 6. Brush the teeth in these sections. Start with the outer surface. Move the brush head across the teeth in small, oval, back-and-forth motions.
⚠ CAUTION: Do not maintain jetting at the same interdental gap for an extended period.
 7. In Auto-jet mode, the camera will detect interdental gaps and trigger a jet. You can manually trigger a jet by pressing and holding the Jet button.
 8. Press the Power button to pause the toothbrush before spitting out the excess fluid. Return the brush head in the mouth and press the Power button again to continue brushing.
⚠ CAUTION: Remove the brush head from the mouth before spitting out excess fluid.
 9. Press the Power button at the end of the brushing session to end operation of the device. It will remain in standby, ready for your next use. You can activate Sleep mode by holding the Power button for 2 seconds.
-

Charging off the dock 5

1. Before charging, ensure there is no moisture, dust or dirt on the charging contacts of the device and the charging cable.
 2. Plug the charging cable into a listed/certified USB-C power adapter and connect the cable to the charging contacts of the toothbrush. Place the device in an upright position.
 3. Check the battery indicator to confirm the device is charging.
 4. To remove the charging cable from the handle, hold the connector and gently pull the brush handle away.
-

Refill dock

The Refill dock provides a convenient way to store and charge your toothbrush and refill the tank. You can choose to use the device without the dock, by using the charging cable provided and refilling the tank manually.

Setup 6

1. Uncoil the USB-C cable from the storage underneath the dock base and connect to a listed/certified USB-C power adapter (5V DC, 500mA–3A), ensuring the charging cable is dry and the dock is upright.
2. Align and attach the base plate to the base of the dock.

The device can be used with or without the Refill dock, which provides convenient charging and tank filling.

Refilling the Dock reservoir **7**

1. Remove the Dock reservoir's cap by pulling up.
2. Fill with Dyson mouth rinse where available, or an alternative low-foaming mouth rinse.
 **CAUTION:** Use only mouth rinse or water.
Do not use hot liquid.
3. Push the Dock reservoir's cap on until it is tightly secured.

Charging and refilling the toothbrush on the Refill dock **8**

1. Before charging, ensure there is no moisture, dust or dirt on the charging contacts of your toothbrush and the charging pins on the Refill dock. Wipe off any moisture, dust or dirt.
2. Attach the brush head and toothbrush tank to the toothbrush handle.
3. Place the toothbrush into the dock base, so that the buttons on your toothbrush face away from the dock.
4. Push the device down into the dock base, insert the tank refilling inlet holes into the nozzle on the dock, and release.
5. Your toothbrush will charge if needed.
6. To refill your toothbrush's tank, press and hold the Dock refill button on the top of the dock. Overfill will be channelled back into the dock reservoir.
7. To remove your toothbrush and stop charging, push the device downwards, against the spring pad to remove it from the dock.

Storage of your device

To store your toothbrush within the Refill dock, empty the refill dock, connect it to the dock as above, and slide the Travel case over the dock until it clicks in place.

You can connect your device to the MyDyson™ app using Wi-Fi and Bluetooth® to access step-by-step support, live cleaning footage, guided cleaning and coverage mapping.

General cleaning instructions **9**

After each use

1. Remove the brush head and rinse under running water.
2. Remove the tank without twisting.
3. Rinse the tank.
4. Clean external surfaces of the device by rinsing.
5. Press and hold the Jet button until no liquid jets out.
6. Dry external surfaces of the device.
7. Allow brush head to dry fully between uses.

Weekly cleaning instructions

It is recommended to clean weekly as per the following steps. If the device is exposed to mould or develops an odour, you should replace the part.

Do not clean the tank and internal parts using harsh chemicals. Ensure the device is turned off and not being charged.

1. Remove the brush head from the handle.
2. Remove the tank by firmly pulling away from the handle, without twisting.
3. Rinse the brush head, nozzle tip and entire handle thoroughly under running water, especially the inlet filter, the brush head connection and the charging contacts.
4. Face the toothbrush into the sink and press and hold the Jet button until no liquid comes out of the nozzle tip.
5. Wipe all external surfaces of the toothbrush with a damp cloth and leave them to air dry.

To clean the tank **10**

1. Remove the tank from the device.
2. Pour hot, boiled water into a suitable container.
3. Use a heatproof utensil to carefully submerge the tank in the water.
4. Ensure the tank fills by submerging for at least 30 seconds.
5. Allow the water to cool before removing the tank.
6. Gently shake the excess water from the tank and leave it to air dry.
7. Reattach the tank by aligning the inlet and firmly pressing the tank into place.

 **CAUTION:** Be careful not to scald yourself when removing the tank from the water.

 **CAUTION:** Only boil water. Do not boil any other liquid.

Removing hard water deposits from the toothbrush **11**

We recommend using Dyson mouth rinse. However, where these are unavailable or water is used in the device, hard water deposits may build up, depending on the water hardness. This may affect the performance of the device. If in a hard water area, follow the instructions below monthly.

1. Remove the tank and brush head from the device.
2. Prepare a container with white vinegar (with 5% acidity) and soak the nozzle tip in it for 5 minutes.
3. Rinse the nozzle tip under running water.
4. Wipe the nozzle tip dry.
5. Fill the tank with white vinegar (with 5% acidity).
6. Attach the tank to the toothbrush handle.
7. Face the toothbrush into the sink, then press and hold the Jet button repeatedly until tank is empty and the device no longer jets liquid.
8. Remove the tank and empty the leftover solution.
9. Repeat steps 5 to 8 with drinkable water instead of vinegar. Do this at least 2 times and until there is no more vinegar residue.
10. Wipe external surfaces of the tank and leave it to air dry.
11. Leave the toothbrush to air dry.

Cleaning the dock 12

It is recommended to clean the dock weekly.

1. Remove the toothbrush from the dock.
2. Remove the Dock reservoir cap.
3. Empty the Dock reservoir.
4. Place a container at the dock base, where the toothbrush would be.
5. Double press the Refill button to eject all the water from the reservoir.
6. Disconnect the dock from the power supply.
7. Use mild soap with warm water and a soft bottle brush to gently clean the reservoir.
8. Rinse the reservoir with running water and wipe external surfaces dry.
9. Clean the charging pins and the dock base with a damp cloth and wipe dry.
10. Wipe off any dirt or buildup at the nozzle.

 **CAUTION:** Clean only the accessible areas.

Do not attempt to disassemble the dock

 **CAUTION:** Do not immerse the Refill dock in water.

 **CAUTION:** Do not pour boiling water into the Dock reservoir.

Removing hard water deposits from the Refill dock

1. Fill up to quarter of the reservoir with white vinegar (with 5% acidity).
2. Place a container at the dock base.
3. Double press the dock button to eject all the vinegar from the dock reservoir.
4. Repeat steps 1 to 3 with drinkable water instead of vinegar. Do this at least 2 times and until there is no more vinegar residue.
5. Leave the dock to air dry.

Environment requirements

Operating conditions: Temperature: 5 °C to 35 °C. Relative Humidity: 15%rH to 90%rH
Atmospheric Pressure: 700hPa to 1060hPa

Transport & Storage conditions:

Temperature: -25 °C to 60 °C. Relative Humidity: 10%rH to 95%rH

Atmospheric Pressure: 500hPa to 1060hPa

Compliance information

The device operates within the following frequency band(s): Bluetooth low energy and 802.11b/g/n. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the device.

Hereby, Dyson Operations Pte. Ltd. declares that the radio equipment type "OC01" is in compliance with Directive 2014/53/EU. The full texts of the declarations of conformity is available at the following internet address:

www.dyson.co.uk/inside-dyson/terms/compliance

- 2400 MHz - 2483.5 MHz; 100 mW EIRP max
- 13.56MHz; 60 dBuA/m max at 10m
- Device networked standby power: < 2 W

To disable or enable Bluetooth low energy and Wi-Fi connectivity, press and hold the Intensity modes button for 4 seconds when the device is on charge.

Limited 2-year guarantee

Terms and conditions of the Dyson 2-year limited guarantee are set out below. The name, address details of the Dyson guarantor for your country are set out in this document – please refer to the table on the back page for full details.

The Dyson guarantee provides benefits which are additional to and do not affect your legal rights as a consumer. The guarantee will apply whether you purchased your device directly from Dyson or from a third party.

What is covered by the guarantee

- The repair or replacement of your Dyson device (at Dyson's discretion) if it is found to be defective due to faulty materials, workmanship or function within two years of purchase or delivery.
- Where this device is sold outside the EU: this guarantee remains valid if the device is used abroad but may only be exercised in the country in which the device is purchased. Please contact Dyson for further information.
- Where this device is sold within the EU: this guarantee remains valid if the device is used abroad. This guarantee may be exercised (i) in the country in which the device was purchased or, (ii) in any of the EU member States for which a Dyson entity is listed as guarantor on the table at the end of this leaflet provided that the same model as this device is sold in the relevant country. Please contact Dyson for further information.

What is not covered by the guarantee

The guarantee does not cover normal wear and tear of the product and any consumable items or components (such as brush heads), or the normal expected reduction in battery discharge time as the battery ages.

In addition, the repair or replacement of a device is not covered where a defect is the result of:

- Damage caused by not carrying out the recommended device maintenance.
- Accidental damage, faults caused by negligent use or care, misuse, neglect, carelessness or operation or handling of the device which is not in accordance with the Dyson User manual.
- Use of the device for professional use or any activity outside of normal domestic use.
- Use of the device for anything other than the use it is designed for.
- Use of parts not assembled or installed in accordance with the instructions of Dyson.
- Use of parts and accessories which are not genuine Dyson components.
- Repairs or alterations carried out by parties other than Dyson or its authorised agents.

If you are in any doubt as to what is covered by your guarantee, please contact the Dyson Helpline.

Summary of cover

- The guarantee becomes effective from the date of purchase (or the date of delivery if this is later).
- Except for devices purchased directly from Dyson, you must provide proof of (both the original and any subsequent) delivery/purchase before any work can be carried out on your Dyson device. Without this proof, any work carried out will be chargeable. Keep your receipt or delivery note.
- All work will be carried out by Dyson or its authorised agents.
- Any parts which are replaced by Dyson will become the property of Dyson.
- The repair or replacement of your Dyson device under guarantee will not extend the period of guarantee unless this is required by local law in the country of purchase.

Important data protection information

You will need to provide us with basic contact information when you register your Dyson device or via the MyDyson™ app. In addition, data including camera data is processed when you use the device.

When registering your Dyson device:

- You will need to provide us with basic contact information to register your device and enable us to support your guarantee.

When registering via the MyDyson™ app:

- You will need to provide us with basic contact information to register the MyDyson™ app. This enables us to securely link your device to your instance of the app.
- When you register, you will have the opportunity to choose whether you would like to receive communications from us. If you opt-in to communications from Dyson, we will send you details of special offers and news of our latest innovations.

When using the device:

Camera: When Auto-jet is on, and the indicator is green, the camera activates to automatically detect gaps and trigger the Conical jet. The camera is turned off at all other times. All images from the camera are processed directly on the device and are not stored. Live images never leave the device, unless you choose to enable live cleaning footage in the app. Even then, images are sent straight from the device to the MyDyson™ app and are deleted immediately after viewing. Dyson and third parties cannot access these images at any time.

Control and explicit consent for oral health: Information relating to your oral health will only be processed with your explicit consent, which you can withdraw at any time via the MyDyson™ app. Certain other data may be processed on other lawful bases to support the safe and effective operation of your device, provide features, or improve your experience. More information can be found in our Privacy Policy.

We never sell your information to third parties. Full details on how we handle your data and your rights are defined by our privacy policies which are available on our website at [privacy.dyson.com](https://www.dyson.com/privacy).

Contact

Manufacturer Name: DYSON OPERATIONS PTE. LTD.
Manufacturer Address: 3 Sentosa Gateway, Singapore 098544
Date of Issue: 18/06/2026
Product Model Reference: OC01

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Notkun 4

1. Facing the brush head into a sink, press and hold the Jet button until all air is expelled and the nozzle begins to jet fluid.
 2. If you are using Brush and Auto-jet mode, or Brush mode, apply toothpaste. Ensure the camera lens is not covered. If you are using Brush and Auto-jet mode, or Brush mode, apply toothpaste. Use the Dyson toothpaste where available, or a non-foaming alternative. Ensure the camera lens is not covered. Place the toothbrush head flat against your teeth. Hold horizontally to follow the jawline.
 3. Place the toothbrush head flat against your teeth. Hold horizontally to follow the jawline.
 4. To begin operation, press the Power button.
 5. Use light pressure on the teeth and gums and let the toothbrush do the work.
 6. Guide the toothbrush slowly across the teeth, moving from tooth to tooth with a gentle oval motion. Do not scrub.
 7. Brush the teeth in these sections. Start with the outer surface. Move the brush head across the teeth in small, oval, back-and-forth motions.
- VARÚÐ:** Ekki halda bunu á sama tannmillibili í langan tíma.
8. Í sjálfvirkri bunustillingu (Auto-jet mode) nemur myndavélin bíl á milli tanna og ræsir bunu. Þú getur ræst bunu með því að ýta á bunuhnappinn og halda honum inni.

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