MOVA Fresh 20 Electric Toothbrush

User Manual

The illustrations in this manual are for reference only.
The actual appearance of the product may vary from what is shown.

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Safety Tips

Safety Warnings

- Toothbrushes can be used by children and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use the appliance in a safe way and understand the hazards involved. Cleaning and user maintenance shall not be made by children without supervision
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children shall not play with the appliance. Cleaning and use maintenance shall not be made by children without supervision.
- The appliance is not intended for use by individuals (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning the use of the appliance by a person responsible for their safety.
- Children should be supervised to ensure they do not play with the appliance.
- This appliance contains a non-replaceable battery. If the battery level remains consistently low, please seek professional assistance. For proper disposal, follow the battery removal instructions.
- Do not put the device or battery in fire or expose them to high temperatures to prevent battery leakage or overheating explosion.
- Please charge the product within a temperature range of 0°C to 45°C.
- Use an adapter with 5V⎓ output voltage to charge this product.
- Unauthorized user modification or repair of this product is strictly prohibited. If the product or power cord is damaged or battery leakage is detected, discontinue use immediately. To avoid danger, please let professionals of the manufacturer handle it.
- Do not use paint thinner, benzene, alcohol and other chemicals to clean this product.
- This product is a Class III appliance. The maximum measured output power of the power supply (power bank) is less than 15W. When connecting the power supply (power bank) to the charging cable for charging the appliance, its rated voltage should be 5V⎓, rated output current 1A, and it must simultaneously meet the requirements of Safety Extra-Low Voltage (SELV).

Precautions for Electric Toothbrush

- When using this product for the first time, you may experience slight bleeding or tooth discomfort, as your gums have not yet adapted to the efficient cleaning. You can start with the Sensitive Mode or customize a gentler brushing mode. This phenomenon typically resolves within 1-2 weeks.
- If bleeding persists for more than 2 weeks, please stop using this product and consult your dentist.

Safety Tips

Precautions for Electric Toothbrush

- Users with severe periodontal disease or those undergoing dental repairs, orthodontics or implants should use this product under the guidance of a dentist.
- This product complies with safety standards for electromagnetic devices. For any questions related to pacemakers or other implanted devices, please consult your physician or the manufacturer of the implanted device before use.
- This product can only be used with the matching brush head.
- Using other brush heads may damage the tooth surface and gums, or affect the lifespan of the electric toothbrush. It is recommended to replace the brush head every 3 months to ensure the cleaning ability of the toothbrush and the cleanliness of the brush head.

Maintenance

- Clean the brush head thoroughly after each use to prevent bacterial growth.
- Replace the brush head regularly; every 3 months is recommended.
- Avoid storing the toothbrush in a damp environment for extended periods to prevent bacteria buildup and battery damage.
- Avoid exposing the toothbrush to extreme heat or cold to protect battery life.
- Do not drop or impact the toothbrush to avoid internal damage.
- Regularly clean the casing and the connection area between the brush head and the body to prevent dirt buildup.
- If not in use for an extended period, clean all components and store the toothbrush in a cool, dry place away from direct sunlight.

Battery Removal

When disposing of this product, the built-in rechargeable battery must be removed and recycled in accordance with local laws and regulations.
Dismantling the battery is not recommended under other circumstances.


1. Ensure the device is powered off and the battery is fully discharged before removing the battery.
2. Use a tool to open the bottom cover of the electric toothbrush and detach it from the main body.
3. Push inward on the latching points on both sides to remove the entire battery assembly.
4. Cut the connection cable between the battery and the motherboard body, take out the battery, and dispose of it at a government-designated location.


Product Introduction

Packing List

Electric toothbrush handle × 1, Whitening brush head × 1, Sensitive care brush head × 1, Toothbrush holder × 1, Type-C charging cable × 1, User manual × 1

Part Description



Function and Usage


Instructions for Use

① Insert the brush head into the handle.
② For first use, long press the power button to exit flight mode. For daily use, pick up the toothbrush to activate.

③ Start/Pause: Short press the power button.
Switch mode: Short press the mode button.
*After 10 seconds of inactivity, the toothbrush will automatically enter sleep mode.

Smart Mode
Adapt to your personal brushing habits through AI-powered learning


Sensitive Mode
Ideal for first-time users and those with sensitive teeth or gums

Cleaning Mode
Suitable for daily use and effective overall cleaning

Whitening Mode
Suitable for users aiming to achieve visibly whiter teeth

Smart Light Strip


- Brushing progress display: During brushing, the light strip gradually fills with color to track your brushing duration.
 *Colors may vary in actual display.
- Example (Based on 120s of brushing time):
 The light strip displays as follows:


Function and Usage

Smart Light Strip

- Smart pressure alert: During brushing, if the brushing pressure is too high, the light strip turns red; if the pressure is too light, it turns yellow.
 *Colors may vary in actual display.


Excessive Pressure


Insufficient Pressure

Smart Display

① During brushing, the screen displays real-time brushing status and reminders.
*The actual content display may vary.


② After brushing, the screen will show the corresponding brushing report to help users immediately learn the brushing status and gradually develop the correct brushing habits.

Current brushing report

Yellow areas indicate inadequately cleaned zone. Rebrush those areas promptly


Plaque report

Visible when using a Plaque Visualization Detector

Function and Usage

Function Settings

- Long press the mode button to enter the setting menu, where you can adjust various functions and configure the on/off settings.
 Short press the Power Button: Confirm and enter settings
 Short press the Mode Button: Next

Adjust Brushing Parameters

- In the setting menu, users can adjust the brushing mode and corresponding parameters (vibration intensity, sweeping angle, and sweeping speed) to suit their brushing needs.


Adjust Brushing Duration

- In the setting menu, users can adjust the brushing duration to suit their needs.

Flight Mode

- In the setting menu, users can enter the Flight Mode. In flight mode, the toothbrush cannot be activated by shaking or short pressing buttons.
- The toothbrush will automatically activate flight mode if shaken multiple times without initiating a brushing session.
- To exit flight mode, long press the power button.

Function and Usage

Bone Conduction Voice Activated Switch

- Enable or disable the voice prompt feature in the setting menu.

Language Settings


- Switch between languages in the setting menu.

Brush Head Usage Reminder

- You can view the number of remaining days for the brush head in the setting menu. After replacing the brush head, you can reset the counter.

System Reset


- System reset clears data and unlink the device from your mobile phone.


Function and Usage

Battery Level Alert

- Charge the toothbrush promptly when the battery is low to prevent it from failing to power on.

Battery Low
Low power, but functional

Insufficient Battery
Charging required before use

APP Interaction

- Scan the QR code on the screen guide page or the one below to download the MOVAhome APP. Make sure the toothbrush is powered on, and follow the APP instructions to pair the toothbrush with the APP.
- In the APP, you can not only adjust ten-level brushing parameters across three dimensions based on personal preferences, but also check brushing reports and configure settings.


Scan the QR code to download MOVAhome APP

Charging

- Open the waterproof cover at the bottom. If moisture is present at the charging port, wipe it dry.
- Gently open the waterproof cover at the bottom of the product. Do not pull it axially with force.
- Use the supporting Type-C charging cable and a power adapter or power bank that meets national standards. Do not use poor-quality chargers.
- Do not use the toothbrush while charging.
- Do not immerse the toothbrush in water or any other liquid while charging.

Basic Parameters

Product name	MOVA Fresh 20 Electric Toothbrush		
Model	ATB25B		
Rated input voltage	5V ⎓	Rated input power	5 W
Waterproof rating	IPX8	Battery capacity	2000 mAh

CE

We, XingKuang Innovation Technology (Suzhou) Co., Ltd., hereby declare that this equipment is in compliance with the applicable Directives and European Norms, and amendments.
The full text of the EU declaration of conformity is available at the following internet address: www.mova.tech

This marking indicates that this product should not be disposed with other household wastes throughout the EU.
To prevent possible harm to the environment or human health from uncontrolled waste disposal, recycle it responsibly to promote the sustainable reuse of material resources.
To return your used device, please use the return and collection systems or contact the retailer where the product was purchased. They can take this product for environmentally safe recycling.
This appliance is not intended for use by persons (including children) with reduced physical,sensory or mental capabilities, or lack of experience and knowledge, unless they have beengiven supervision or instruction concerning use of the appliance by a person responsible fortheir safety.

Children should be supervised to ensure that they do not play with the appliance.

This appliance contains batteries that are non-replaceable.

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Scan the QR code to download the MOVAhome APP Connect to MOVA electric toothbrush

Scan the QR code to view the product tutorial and electronic manual

22
PAP

User manual version number: V1.3
Service website: www.mova.tech
Manufacturer: XingKuang Innovation Technology (Suzhou) Co., Ltd.
Address:12th Floor, Building 2, Wisdom Valley Park, Suzhou (Taihu) Software Industrial Park, No. 1463 Wuzhong Avenue, Yuexi Street, Wuzhong Economic Development Zone, Suzhou, Jiangsu Province, P.R. China.
Made in China